

TUFT1 antibody

Product Information

Catalog No.:	FNab09106
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in the mineralization and structural organization of enamel.

Immunogen information

Immunogen:	tuftelin 1
Synonyms:	Tuftelin TUFT1
Observed MW:	50 kDa
Uniprot ID :	Q9NNX1

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:5000; IF: 1:20-1:200
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab09106(TUFT1 antibody) at dilution of 1:1000