

TULP4 antibody

Product Information

Catalog No.:	FNab09110
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be a substrate-recognition component of a SCF-like ECS(Elongin-Cullin-SOCS-box protein) E3 ubiquitin ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.

Immunogen information


Immunogen:	tubby like protein 4
Synonyms:	Tubby-related protein 4 Tubby superfamily protein Tubby-like protein 4 TULP4 KIAA1397 TUBL4 TUSP
Observed MW:	169 kDa
Uniprot ID :	Q9NRJ4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab09110(TULP4 antibody) at dilution of 1:100



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab09110(TULP4 antibody) at dilution of 1:800