

## TWNK antibody

### Product Information

Catalog No.:	FNab09115
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	chromosome 10 open reading frame 2
Synonyms:	Twinkle mtDNA helicase Progressive external ophthalmoplegia 1 protein T7 gp4-like protein with intramitochondrial nucleoid localization T7-like mitochondrial DNA helicase Twinkle protein, mitochondrial TWNK C10orf2 PEO1
Observed MW:	60 kDa
Uniprot ID :	Q96RR1

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using FNab09115(PEO1 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09115(PEO1 antibody) at dilution of 1:300