

TXNDC11 antibody

Product Information

Catalog No.:	FNab09122
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May act as a redox regulator involved in DUOX proteins folding. The interaction with DUOX1 and DUOX2 suggest that it belongs to a multiprotein complex constituting the thyroid H₂O₂ generating system. It is however not sufficient to assist DUOX1 and DUOX2 in H₂O₂ generation.

Immunogen information

Immunogen:	thioredoxin domain containing 11
Synonyms:	Thioredoxin domain-containing protein 11 EF-hand-binding protein 1 TXNDC11 EFP1
Observed MW:	30 kDa
Uniprot ID :	Q6PKC3

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:5000
Image:	



DU 145 cells were subjected to SDS PAGE followed by western blot with FNab09122(TXNDC11 Antibody) at dilution of 1:1000