

U2AF2 antibody

Product Information

Catalog No.:	FNab09139
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	U2 small nuclear RNA auxiliary factor 2
Synonyms:	Splicing factor U2AF 65 kDa subunit U2 auxiliary factor 65 kDa subunit (hU2AF(65), hU2AF65) U2 snRNP auxiliary factor large subunit U2AF2 U2AF65
Observed MW:	62 kDa
Uniprot ID :	P26368

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab09139(U2AF2 antibody) at dilution of 1:50



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab09139(U2AF2 antibody) at dilution of 1:1000