

UBAC2 antibody

Product Information

Catalog No.:	FNab09148
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Restricts trafficking of FAF2 from the endoplasmic reticulum to lipid droplets.

Immunogen information

Immunogen:	UBA domain containing 2
Synonyms:	Ubiquitin-associated domain-containing protein 2 (UBA domain-containing protein 2) Phosphoglycerate dehydrogenase-like protein 1 UBAC2 PHGDHL1
Observed MW:	35 kDa
Uniprot ID :	Q8NBM4

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:2000; IF: 1:20-1:200
Image:	



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab09148(UBAC2 Antibody) at dilution of 1:600