

## UBE2D4 antibody

### Product Information

Catalog No.:	FNab09169
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


This antibody can recognize all the four members of UBE2D due to the high homology.

### Immunogen information


Immunogen:	ubiquitin-conjugating enzyme E2 D4(putative)
Synonyms:	Ubiquitin-conjugating enzyme E2 D4 E2 ubiquitin-conjugating enzyme D4 HBUCE1 Ubiquitin carrier protein D4 Ubiquitin-protein ligase D4 UBE2D4 UBCH5D
Observed MW:	17 kDa
Uniprot ID :	Q9Y2X8

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:3000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab09169(UBE2D4 antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab09169(UBE2D4 antibody) at dilution of 1:500