

## UBE2U antibody

### Product Information

Catalog No.:	FNab09186
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)


### Background

#### Immunogen information


Immunogen:	ubiquitin-conjugating enzyme E2U (putative)
Synonyms:	Ubiquitin-conjugating enzyme E2 U E2 ubiquitin-conjugating enzyme U Ubiquitin carrier protein U Ubiquitin-protein ligase U UBE2U
Observed MW:	42 kDa
Uniprot ID :	Q5V VX9

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:100
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNab09186(UBE2U antibody) at dilution of 1:50



human heart tissue were subjected to SDS PAGE followed by western blot with FNab09186(UBE2U antibody) at dilution of 1:1000