

## UBE2Z antibody

### Product Information

Catalog No.:	FNab09188
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background


This gene encodes an enzyme which ubiquitinates proteins which participate in signaling pathways and apoptosis.

### Immunogen information


Immunogen:	ubiquitin-conjugating enzyme E2Z
Synonyms:	Ubiquitin-conjugating enzyme E2 Z E2 ubiquitin-conjugating enzyme Z Uba6-specific E2 conjugating enzyme 1 (Use1) Ubiquitin carrier protein Z Ubiquitin-protein ligase Z UBE2Z
Observed MW:	25 kDa, 39 kDa
Uniprot ID :	Q9H832

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:50 - 1:100
Image:	



Immunohistochemistry of paraffin-embedded mouse kidney using FNab09188(UBE2Z antibody) at dilution of 1:100



human testis tissue were subjected to SDS PAGE followed by western blot with FNab09188(UBE2Z antibody) at dilution of 1:1000