

UBTF antibody

Product Information

Catalog No.:	FNab09194
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	upstream binding transcription factor, RNA polymerase I
Synonyms:	Nucleolar transcription factor 1 Autoantigen NOR-90 Upstream-binding factor 1 (UBF-1) UBTF UBF UBF1
Observed MW:	89 kDa
Uniprot ID :	P17480

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:2000; IF: 1:50-1:500
Image:	



HT-1080 cells were subjected to SDS PAGE followed by western blot with FNab09194(UBF1 Antibody) at dilution of 1:600