

UBL3 antibody

Product Information

Catalog No.:	FNab09198
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	ubiquitin-like 3
Synonyms:	Ubiquitin-like protein 3 Membrane-anchored ubiquitin-fold protein (HsMUB, MUB) Protein HCG-1 UBL3 PNSC1
Observed MW:	13 kDa
Uniprot ID :	O95164

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:5000
Image:	



human brain tissue were subjected to SDS PAGE followed by western blot with FNab09198(UBL3 antibody) at dilution of 1:1500