

## UBXN4 antibody

### Product Information

Catalog No.:	FNab09213
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	UBX domain protein 4
Synonyms:	UBX domain-containing protein 4 Erasin UBX domain-containing protein 2 UBXN4 KIAA0242 UBXD2 UBXDC1
Observed MW:	60 kDa
Uniprot ID :	Q92575

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:500-1:3000; IF:1:10-1:100; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart using FNab09213(UBXN4 antibody) at dilution of 1:100



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab09213(UBXN4 antibody) at dilution of 1:500