

UCHL3 antibody

Product Information

Catalog No.:	FNab09219
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

The protein encoded by this gene is a member of the deubiquitinating enzyme family. Members of this family are proteases that catalyze the removal of ubiquitin from polypeptides and are divided into five classes, depending on the mechanism of catalysis. This protein may hydrolyze the ubiquitinyl-N-epsilon amide bond of ubiquitinated proteins to regenerate ubiquitin for another catalytic cycle. Alternative splicing results in multiple transcript variants that encode different protein isoforms.


Immunogen information

Immunogen:	ubiquitin carboxyl-terminal esterase L3 (ubiquitin thioesterase)
Synonyms:	Ubiquitin carboxyl-terminal hydrolase isozyme L3 (UCH-L3) Ubiquitin thioesterase L3 UCHL3
Observed MW:	26 kDa
Uniprot ID :	P15374


Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:100

Image:



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab09219(UCHL3 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09219(UCHL3 antibody) at dilution of 1:1000