

UCKL1 antibody

Product Information

Catalog No.:	FNab09224
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May contribute to UTP accumulation needed for blast transformation and proliferation.

Immunogen information


Immunogen:	uridine-cytidine kinase 1-like 1
Synonyms:	Uridine-cytidine kinase-like 1 UCKL1 URKL1
Observed MW:	46 kDa, 61 kDa
Uniprot ID :	Q9NWZ5

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200;IP: 1:200-1:1000
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab09224(UCKL1 antibody) at dilution of 1:100



IP Result of anti-UCKL1 (IP:FNab09224, 4ug; Detection:FNab09224 1:500) with K-562 cells lysate 3200ug.



K-562 cells were subjected to SDS PAGE followed by western blot with FNab09224(UCKL1 antibody) at dilution of 1:300