

UGGT2 antibody

Product Information

Catalog No.:	FNab09237
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	UDP-glucose ceramide glucosyltransferase-like 2
Synonyms:	UDP-glucose:glycoprotein glucosyltransferase 2 (UGT2, hUGT2) UDP--Glc:glycoprotein glucosyltransferase 2 UDP-glucose ceramide glucosyltransferase-like 1 UGGT2 UGCGL2 UGT2
Observed MW:	175 kDa
Uniprot ID :	Q9NYU1

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:5000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab09237(UGCGL2 antibody) at dilution of 1:300