

UGT3A2 antibody

Product Information

Catalog No.:	FNab09242
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


UDP-glucuronosyltransferases catalyze phase II biotransformation reactions in which lipophilic substrates are conjugated with glucuronic acid to increase water solubility and enhance excretion. They are of major importance in the conjugation and subsequent elimination of potentially toxic xenobiotics and endogenous compounds(By similarity).

Immunogen information

Immunogen:	UDP glycosyltransferase 3 family, polypeptide A2
Synonyms:	UDP-glucuronosyltransferase 3A2 (UDPGT 3A2) UGT3A2
Observed MW:	55 kDa
Uniprot ID :	Q3SY77

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:5000
Image:	



human testis tissue were subjected to SDS PAGE followed by western blot with FNab09242(UGT3A2 antibody) at dilution of 1:500