

BLTP3B antibody

Product Information

Catalog No.:	FNab09246
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	UHRF1 binding protein 1-like
Synonyms:	Bridge-like lipid transfer protein family member 3B Syntaxin-6 Habc-interacting protein of 164 kDa UHRF1-binding protein 1-like BLTP3B KIAA0701 SHIP164 UHRF1BP1L
Observed MW:	164 kDa
Uniprot ID :	A0JNW5

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IP, WB
Recommended dilution:	WB: 1:500-1:5000; IP: 1:500-1:5000
Image:	



IP Result of anti-UHRF1BP1L (IP:FNab09246, 5ug; Detection:FNab09246 1:1000) with mouse brain tissue lysate 5000ug.



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab09246(UHRF1BP1L Antibody) at dilution of 1:1500