

## ULBP2 antibody

### Product Information

Catalog No.:	FNab09250
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background

Ligand for the KLRK1/NKG2D receptor, together with at least ULBP1 and ULBP3. ULBPs activate multiple signaling pathways in primary NK cells, resulting in the production of cytokines and chemokines. Binding of ULBPs ligands to KLRK1/NKG2D induces calcium mobilization and activation of the JAK2, STAT5, ERK and PI3K kinase/Akt signal transduction pathway. In CMV infected cells, interacts with soluble CMV glycoprotein UL16. The interaction with UL16 blocked the interaction with the KLRK1/NKG2D receptor, providing a mechanism by which CMV infected cells might escape the immune system. UL16 also causes ULBP2 to be retained in the ER and cis-Golgi apparatus so that it does not reach the cell surface.

### Immunogen information


Immunogen:	UL16 binding protein 2
Synonyms:	UL16-binding protein 2 ALCAN-alpha NKG2D ligand 2 (N2DL-2, NKG2DL2) Retinoic acid early transcript 1H ULBP2 N2DL2 RAET1H
Observed MW:	30-35 kDa
Uniprot ID :	Q9BZM5

### Application

Reactivity:	Human
Tested Application:	ELISA, WB

Recommended dilution: WB: 1:500-1:5000

Image:



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab09250(ULBP2 antibody) at dilution of 1:400