

UNC5CL antibody

Product Information

Catalog No.:	FNab09262
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Inhibits NF-kappa-B-dependent transcription by impairing NF-kappa-B binding to its targets.

Immunogen information

Immunogen:	unc-5 homolog C(C. elegans)-like
Synonyms:	UNC5C-like protein Protein unc-5 homolog C-like ZU5 and death domain-containing protein UNC5CL ZUD
Observed MW:	70 kDa
Uniprot ID :	Q8IV45

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:1000
Image:	



mouse colon tissue were subjected to SDS PAGE followed by western blot with FNab09262(UNC5CL Antibody) at dilution of 1:300