

## UQCRB antibody

### Product Information

Catalog No.:	FNab09279
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


This is a component of the ubiquinol-cytochrome c reductase complex(complex III or cytochrome b-c1 complex), which is part of the mitochondrial respiratory chain. This component is involved in redox-linked proton pumping.

### Immunogen information


Immunogen:	ubiquinol-cytochrome c reductase binding protein
Synonyms:	Cytochrome b-c1 complex subunit 7 Complex III subunit 7 Complex III subunit VII QP-C Ubiquinol-cytochrome c reductase complex 14 kDa protein UQCRB UQBP
Observed MW:	48 kDa
Uniprot ID :	P14927

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:1000; IHC:1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab09279(UQCRB antibody) at dilution of 1:100



Various lysates were subjected to SDS PAGE followed by western blot with FNab09279(UQCRB antibody) at dilution of 1:400