

## USF2 antibody

### Product Information

Catalog No.:	FNab09298
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Transcription factor that binds to a symmetrical DNA sequence(E-boxes)(5'-CACGTG-3') that is found in a variety of viral and cellular promoters.

### Immunogen information

Immunogen:	upstream transcription factor 2, c-fos interacting
Synonyms:	Upstream stimulatory factor 2 Class B basic helix-loop-helix protein 12 (bHLHb12) FOS-interacting protein (FIP) Major late transcription factor 2 Upstream transcription factor 2 USF2 BHLHB12
Observed MW:	44 kDa
Uniprot ID :	Q15853

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:5000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab09298(USF2 antibody) at dilution of 1:1000