

USHBP1 antibody

Product Information

Catalog No.:	FNab09300
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	Usher syndrome 1C binding protein 1
Synonyms:	Harmonin-binding protein USHBP1 AIE-75-binding protein Mutated in colon cancer protein 2 (MCC-2) Usher syndrome type-1C protein-binding protein 1 (USH1C-binding protein 1) USHBP1 AIEBP MCC2
Observed MW:	90 kDa
Uniprot ID :	Q8N6Y0

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human heart tissue slide using FNab09300(USHBP1 Antibody) at dilution of 1:50



A549 cells were subjected to SDS PAGE followed by western blot with FNab09300(USHBP1 Antibody) at dilution of 1:300