

anti- USP18 antibody

Product Information

Catalog No.:	FNab09313
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Can efficiently cleave only ISG15 fusions including native ISG15 conjugates linked via isopeptide bonds. Necessary to maintain a critical cellular balance of ISG15-conjugated proteins in both healthy and stressed organisms. Isoform 2: Has enzymatic activity similar to isoform 1 and interferes with type I interferon signaling. Major deISGylation enzyme for nuclear proteins(PubMed:22170061).

Immunogen information

Immunogen:	ubiquitin specific peptidase 18
Synonyms:	43 kDa ISG15 specific protease, hUBP43, ISG43, Ubl thiolesterase 18, UBP43, USP18
Observed MW:	70 kDa, 43 kDa
Uniprot ID :	Q9UMW8

Application

Reactivity:	Human
Tested Application:	ELISA, WB

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com

Recommended dilution: WB: 1:500-1:1000

Image:



WB result of USP18 antibody(FNab09313) with MCF-7 cells .

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com