

UTP3 antibody

Product Information

Catalog No.:	FNab09350
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Essential for gene silencing: has a role in the structure of silenced chromatin. Plays a role in the developing brain(By similarity).

Immunogen information

Immunogen:	UTP3, small subunit(SSU) processome component, homolog
Synonyms:	Something about silencing protein 10 Charged amino acid-rich leucine zipper 1 (CRL1) Disrupter of silencing SAS10 UTP3 homolog UTP3 CRLZ1 SAS10
Observed MW:	30 kDa, 70 kDa
Uniprot ID :	Q9NQZ2

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:1000-1:10000; IF: 1:20-1:200
Image:	



A549 cells were subjected to SDS PAGE followed by western blot with FNab09350(UTP3 antibody) at dilution of 1:1000