

UTP6 antibody

Product Information

Catalog No.:	FNab09351
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in nucleolar processing of pre-18S ribosomal RNA.

Immunogen information

Immunogen:	UTP6, small subunit(SSU) processome component, homolog(yeast)
Synonyms:	U3 small nucleolar RNA-associated protein 6 homolog Hepatocellular carcinoma-associated antigen 66 Multiple hat domains protein UTP6 C17orf40 HCA66 MHAT
Observed MW:	78 kDa, 36 kDa
Uniprot ID :	Q9NYH9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200-1:1000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09351(UTP6 Antibody) at dilution of 1:300