

## VAMP2 antibody

### Product Information

Catalog No.:	FNab09359
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Involved in the targeting and/or fusion of transport vesicles to their target membrane. Modulates the gating characteristics of the delayed rectifier voltage-dependent potassium channel KCNB1.

### Immunogen information


Immunogen:	vesicle-associated membrane protein 2(synaptobrevin 2)
Synonyms:	Vesicle-associated membrane protein 2 (VAMP-2) Synaptobrevin-2 VAMP2 SYB2
Observed MW:	18 kDa
Uniprot ID :	P63027

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human skeletal muscle using FNab09359(VAMP2 antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab09359(VAMP2 Antibody) at dilution of 1:600