

VAMP5 antibody

Product Information

Catalog No.:	FNab09363
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May participate in trafficking events that are associated with myogenesis, such as myoblast fusion and/or GLUT4 trafficking.

Immunogen information


Immunogen:	vesicle-associated membrane protein 5(myobrevin)
Synonyms:	Vesicle-associated membrane protein 5 (VAMP-5) Myobrevin VAMP5
Observed MW:	13 kDa
Uniprot ID :	O95183

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNab09363(VAMP5 antibody) at dilution of 1:50



human lung tissue were subjected to SDS PAGE followed by western blot with FNab09363(VAMP5 antibody) at dilution of 1:400