

VAPB antibody

Product Information

Catalog No.:	FNab09371
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	9C9
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Participates in the endoplasmic reticulum unfolded protein response(UPR) by inducing ERN1/IRE1 activity. Involved in cellular calcium homeostasis regulation.

Immunogen information


Immunogen:	VAMP(vesicle-associated membrane protein)-associated protein B and C
Synonyms:	Vesicle-associated membrane protein-associated protein B/C (VAMP-B/VAMP-C, VAMP-associated protein B/C, VAP-B/VAP-C) VAPB
Observed MW:	27 kDa
Uniprot ID :	O95292

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200;IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab09371(VAPB Antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab09371(VAPB Antibody) at dilution of 1:1000