

## VASH2 antibody

### Product Information

Catalog No.:	FNab09374
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Angiogenesis inhibitor. Inhibits network formation by endothelial cells.

### Immunogen information

Immunogen:	vasohibin 2
Synonyms:	Tubuliny1-Tyr carboxypeptidase 2 Vasohibin-2 Vasohibin-like protein VASH2 VASHL
Observed MW:	40 kDa
Uniprot ID :	Q86V25

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:5000; IF: 1:10-1:100
Image:	



A549 cells were subjected to SDS PAGE followed by western blot with FNab09374(VASH2 antibody) at dilution of 1:600