

VGLL1 antibody

Product Information

Catalog No.:	FNab09398
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May act as a specific coactivator for the mammalian TEFs.

Immunogen information


Immunogen:	vestigial like 1(Drosophila)
Synonyms:	Transcription cofactor vestigial-like protein 1 (Vgl-1) Protein TONDU VGLL1 TDU
Observed MW:	29 kDa
Uniprot ID :	Q99990

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human prostate cancer using FNab09398(VGLL1 antibody) at dilution of 1:50



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab09398(VGLL1 antibody) at dilution of 1:1000