

VPS72 antibody

Product Information

Catalog No.:	FNab09450
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Could be a DNA-binding transcriptional regulator. May be involved in chromatin modification and remodeling.

Immunogen information

Immunogen:	vacuolar protein sorting 72 homolog
Synonyms:	Vacuolar protein sorting-associated protein 72 homolog Protein YL-1 Transcription factor-like 1 VPS72 TCFL1 YL1
Observed MW:	57 kDa
Uniprot ID :	Q15906

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200-1:2000
Image:	



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab09450(VPS72 antibody) at dilution of 1:800