

## VSIG1 antibody

### Product Information

Catalog No.:	FNab09454
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	V-set and immunoglobulin domain containing 1
Synonyms:	V-set and immunoglobulin domain-containing protein 1 Cell surface A33 antigen Glycoprotein A34 VSIG1 GPA34
Observed MW:	70 kDa
Uniprot ID :	Q86XK7

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000;IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab09454(VSIG1 antibody) at dilution of 1:100



human stomach tissue were subjected to SDS PAGE followed by western blot with FNab09454(VSIG1 antibody) at dilution of 1:1500