

VSTM1 antibody

Product Information

Catalog No.:	FNab09456
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Isoform 2: Behaves as a cytokine, promoting IL17A secretion by CD4+ T-cells, and differentiation and activation of IL17 producing helper T-cells(TH17). Isoform 1: Inhibitory immune receptor involved in the regulation of phagocytes.

Immunogen information

Immunogen:	V-set and transmembrane domain containing 1
Synonyms:	V-set and transmembrane domain-containing protein 1 Signal inhibitory receptor on leukocytes-1 (SIRL-1) VSTM1
Observed MW:	40 kDa
Uniprot ID :	Q6UX27

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:1000
Image:	



HL-60 cells were subjected to SDS PAGE followed by western blot with FNab09456(VSTM1 Antibody) at dilution of 1:300