

anti- VSX1 antibody

Product Information

| | |
|---------------|---|
| Catalog No.: | FNab09458 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Background

Binds to the 37-bp core of the locus control region(LCR) of the red/green visual pigment gene cluster. May regulate the activity of the LCR and the cone opsin genes at earlier stages of development.

Immunogen information

| | |
|--------------|--|
| Immunogen: | visual system homeobox 1 |
| Synonyms: | Homeodomain protein RINX, KTCN, PPCD, PPD, RINX, Transcription factor VSX1, visual system homeobox 1, VSX1 |
| Observed MW: | 38 kDa, observed 70 kDa |
| UniprotID : | Q9NZR4 |

Application

| | |
|--------------------------|--------------|
| Reactivity: | Human, Mouse |
| Tested Application: | ELISA, WB |
| Recommended dilution:WB: | 1:200-1:1000 |


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com

Image:



mouse eye tissue were subjected to SDS PAGE followed by western blot with FNab09458(VSX1 Antibody) at dilution of 1:300

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com