

## VSX2 antibody

### Product Information

Catalog No.:	FNab09459
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Plays a significant role in the specification and morphogenesis of the sensory retina. May also participate in the development of the cells of the inner nuclear layer, particularly bipolar cells(By similarity).

### Immunogen information

Immunogen:	visual system homeobox 2
Synonyms:	Visual system homeobox 2 Ceh-10 homeodomain-containing homolog Homeobox protein CHX10 VSX2 CHX10 HOX10
Observed MW:	42 kDa
Uniprot ID :	P58304

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:2000
Image:	



mouse eye tissue were subjected to SDS PAGE followed by western blot with FNab09459(VSX2 Antibody) at dilution of 1:600