

## VWA1 antibody

### Product Information

Catalog No.:	FNab09464
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Promotes matrix assembly.

### Immunogen information


Immunogen:	von Willebrand factor A domain containing 1
Synonyms:	von Willebrand factor A domain-containing protein 1 VWA1
Observed MW:	47 kDa
Uniprot ID :	Q6PCB0

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:200-1:2000;IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab09464(VWA1 antibody) at dilution of 1:50



Recombinant proteins were subjected to SDS PAGE followed by western blot with FNab09464(VWA1 Antibody) at dilution of 1:800