

WASL antibody

Product Information

Catalog No.:	FNab09470
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

Regulates actin polymerization by stimulating the actin-nucleating activity of the Arp2/3 complex. Involved in mitosis and cytokinesis, via its role in the regulation of actin polymerization. Binds to HSF1/HSTF1 and forms a complex on heat shock promoter elements(HSE) that negatively regulates HSP90 expression.

Immunogen information


Immunogen:	Wiskott-Aldrich syndrome-like
Synonyms:	Actin nucleation-promoting factor WASL Neural Wiskott-Aldrich syndrome protein (N-WASP) WASL
Observed MW:	55 kDa
Uniprot ID :	O00401

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200;IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab09470(WASL antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab09470(WASL antibody) at dilution of 1:500