

WBP2 antibody

Product Information

Catalog No.:	FNab09474
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	WW domain binding protein 2
Synonyms:	WW domain-binding protein 2 (WBP-2) WBP2
Observed MW:	28 kDa, 40 kDa
Uniprot ID :	Q969T9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:5000;IHC: 1:20-1:200;IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab09474(WBP2 antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab09474(WBP2 antibody) at dilution of 1:1000