

## WBP2NL antibody

### Product Information

Catalog No.:	FNab09475
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May play a role in meiotic resumption and pronuclear formation, mediated by a WW domain-signaling pathway during fertilization.

### Immunogen information


Immunogen:	WBP2 N-terminal like
Synonyms:	Postacrosomal sheath WW domain-binding protein WW domain-binding protein 2-like WBP2NL PAWP
Observed MW:	28 kDa
Uniprot ID :	Q6ICG8

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:200-1:2000;IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab09475(WBP2NL Antibody) at dilution of 1:50



human testis tissue were subjected to SDS PAGE followed by western blot with FNab09475(WBP2NL Antibody) at dilution of 1:600