

WDR20 antibody

Product Information

Catalog No.:	FNab09486
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	WD repeat domain 20
Synonyms:	WD repeat-containing protein 20 Protein DMR WDR20
Observed MW:	63 kDa
Uniprot ID :	Q8TBZ3

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:5000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09486(WDR20 Antibody) at dilution of 1:1000