

WDR33 antibody

Product Information

Catalog No.:	FNab09490
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Essential for both cleavage and polyadenylation of pre-mRNA 3' ends.

Immunogen information

Immunogen:	WD repeat domain 33
Synonyms:	pre-mRNA 3' end processing protein WDR33 WD repeat-containing protein 33 WD repeat-containing protein of 146 kDa WDR33 WDC146
Observed MW:	150 kDa
Uniprot ID :	Q9C0J8

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:5000; IP: 1:500-1:5000
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab09490(WDR33 Antibody) at dilution of 1:1000