

WDR7 antibody

Product Information

Catalog No.:	FNab09501
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	WD repeat domain 7
Synonyms:	WD repeat-containing protein 7 Rabconnectin-3 beta TGF-beta resistance-associated protein TRAG WDR7 KIAA0541 TRAG
Observed MW:	145 kDa
Uniprot ID :	Q9Y4E6

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis slide using FNab09501(WDR7 Antibody) at dilution of 1:50



rat brain tissue were subjected to SDS PAGE followed by western blot with FNab09501(WDR7 Antibody) at dilution of 1:1000