

WNT5A-B antibody

Product Information

Catalog No.:	FNab09520
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	wingless-type MMTV integration site family, member 5A
Synonyms:	Protein Wnt-5a WNT5A
Observed MW:	42 kDa
Uniprot ID :	P41221

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using FNab09520(WNT5A-B Antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab09520(WNT5A-B Antibody) at dilution of 1:600