

WNT7A antibody

Product Information

Catalog No.:	FNab09522
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. Signaling by Wnt-7a allows sexually dimorphic development of the mullerian ducts(By similarity).

Immunogen information


Immunogen:	wingless-type MMTV integration site family, member 7A
Synonyms:	Protein Wnt-7a WNT7A
Observed MW:	39 kDa
Uniprot ID :	O00755

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:5000; IF: 1:20-1:200; IP: 1:200-1:1000
Image:	



IP Result of anti-WNT7A (IP:FNab09522, 4ug;
Detection:FNab09522 1:300) with mouse kidney
tissue lysate 4000ug.



mouse kidney tissue were subjected to SDS PAGE
followed by western blot with
FNab09522(WNT7A Antibody) at dilution of
1:1000