

XAB2 antibody

Product Information

Catalog No.:	FNab09536
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in transcription-coupled repair(TCR), transcription and pre-mRNA splicing.

Immunogen information


Immunogen:	XPA binding protein 2
Synonyms:	Pre-mRNA-splicing factor SYF1 Protein HCNP XPA-binding protein 2 XAB2 HCNP KIAA1177 SYF1
Observed MW:	100 kDa
Uniprot ID :	Q9HCS7

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:200-1:1000; IP: 1:200-1:1000
Image:	



IP Result of anti-XAB2 (IP:FNab09536, 4ug;
Detection:FNab09536 1:200) with K-562 cells
lysate 3200ug.



K-562 cells were subjected to SDS PAGE
followed by western blot with FNab09536(XAB2
antibody) at dilution of 1:200