

## XIRP2 antibody

### Product Information

Catalog No.:	FNab09542
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background

Protects actin filaments from depolymerization.

### Immunogen information


Immunogen:	xin actin-binding repeat containing 2
Synonyms:	Xin actin-binding repeat-containing protein 2 Beta-xin Cardiomyopathy-associated protein 3 Xeplin XIRP2 CMYA3
Observed MW:	280 kDa
Uniprot ID :	A4UGR9

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human skeletal muscle using FNab09542(XIRP2 antibody) at dilution of 1:50



COLO 320 cells were subjected to SDS PAGE followed by western blot with FNab09542(XIRP2 antibody) at dilution of 1:500