

XK antibody

Product Information

Catalog No.:	FNab09543
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in sodium-dependent transport of neutral amino acids or oligopeptides.

Immunogen information


Immunogen:	Membrane transport protein XK
Synonyms:	Endoplasmic reticulum membrane adapter protein XK Kell complex 37 kDa component Kx antigen Membrane transport protein XK XK-related protein 1 XK XKR1 XRG1
Observed MW:	42 kDa
Uniprot ID :	P51811

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:50-1:400
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer tissue slide using FNab09543(XK Antibody) at dilution of 1:200



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with FNab09543(XK Antibody) at dilution of 1:600