

YIPF4 antibody

Product Information

Catalog No.:	FNab09564
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in the maintenance of the Golgi structure.

Immunogen information


Immunogen:	Yip1 domain family, member 4
Synonyms:	Protein YIPF4 YIP1 family member 4 YIPF4
Observed MW:	28 kDa
Uniprot ID :	Q9BSR8

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:5000; IF: 1:10-1:100; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human endometrial cancer using FNab09564(YIPF4 antibody) at dilution of 1:50



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab09564(YIPF4 antibody) at dilution of 1:500