

## YY1AP1 antibody

### Product Information

Catalog No.:	FNab09579
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Enhances transcription activation by YY1. May play a role in cell cycle regulation.

### Immunogen information

Immunogen:	YY1 associated protein 1
Synonyms:	YY1-associated protein 1 Hepatocellular carcinoma susceptibility protein Hepatocellular carcinoma-associated protein 2 YY1AP1 HCCA2 YY1AP
Observed MW:	60 kDa, 75 kDa
Uniprot ID :	Q9H869

### Application

Reactivity:	Human
Tested Application:	ELISA, IF, WB
Recommended dilution:	WB: 1:200-1:2000; IF: 1:10-1:100
Image:	



MDA-MB-231 cells were subjected to SDS PAGE followed by western blot with FNab09579 (YY1AP1 antibody) at dilution of 1:800