

ZBTB11 antibody

Product Information

Catalog No.:	FNab09588
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in transcriptional regulation.

Immunogen information

Immunogen:	zinc finger and BTB domain containing 11
Synonyms:	Zinc finger and BTB domain-containing protein 11 ZBTB11
Observed MW:	69 kDa, 140 kDa
Uniprot ID :	O95625

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200-1:2000
Image:	



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab09588(ZBTB11 Antibody) at dilution of 1:1000